## **1. Información general:**

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The screenshot displays the Pentacam software interface with data for two eyes: 'Derecho' (Right) and 'Izquierdo' (Left). The interface is divided into several sections: patient information, biometric measurements, and refractive data. A red box highlights the bottom section, which is used for entering optimization data. This section includes fields for 'Fecha de la cirugía' (Surgery date), 'IOL usada' (IOL used), 'Potencia' (Power), and 'Fórmula' (Formula). The 'Base de datos de optimización' (Optimization data base) is also indicated.

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Add New Surgery Details ✕

**Notes:**

- Steep K must be greater than Flat K for Pre-op and Post-op measurements.
- Steep/Flat Meridian must be within 90 +/- 5 degrees for Pre-op and Post-op measurements.
- Potentially incorrect data will be highlighted in orange

Chart Number		Date Of Surgery		Patient Age	
<input type="text"/>		<input type="text" value="Enter Date"/>		<input type="text" value="0"/>	
Eye	<input type="text" value="-Select-"/>	Prior Corneal Surgery	<input type="checkbox"/>	Incision Location (degrees)	<input type="text" value="0-360"/>
Incision Type	<input type="text" value="-Select-"/>	Incision Size (0.01 to 9.99 mm)	<input type="text" value="0.01-9.99"/>		
Pre-Op K1	<input type="text" value="35-50"/>	Pre-Op K2	<input type="text" value="35-50"/>		
Pre-Op A1	<input type="text" value="0-180"/>	Pre-Op A2	<input type="text" value="0-180"/>		
Post-Op K1	<input type="text" value="35-50"/>	Post-Op K2	<input type="text" value="35-50"/>		
Post-Op A1	<input type="text" value="0-180"/>	Post-Op A2	<input type="text" value="0-180"/>		

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- Biometry and Intraocular Lens Power Calculation." *Experimental and Therapeutic Medicine* 21, no. 4 (2021). doi:10.3892/etm.2021.9757.
5. Olsen, Thomas, and Peter Hoffmann. "C constant: New Concept for Ray Tracing-assisted Intraocular Lens Power Calculation." *Journal of Cataract and Refractive Surgery* 40, no. 5 (2014): 764-773. doi:10.1016/j.jcrs.2013.10.037.
  6. Schröder, Sven, Christina Leydolt, Rupert Menapace, Timo Eppig, and Achim Langenbacher. "Determination of Personalized IOL-constants for the Haigis Formula under Consideration of Measurement Precision." *PLoS One* 11, no. 7 (2016). doi:10.1371/journal.pone.0158988.
  7. Aristodemou, Panayiotis, N. E. Knox Cartwright, J. M. Sparrow, and R. L. Johnston. "Intraocular Lens Formula Constant Optimization and Partial Coherence Interferometry Biometry: Refractive Outcomes in 8108 Eyes after Cataract Surgery." *Journal of Cataract and Refractive Surgery* 37, no. 1 (2011): 50-62. doi:10.1016/j.jcrs.2010.07.037.
  8. Charalampidou, Sofia, Lynne Cassidy, Ed Ng, et al. "Effect on Refractive Outcomes after Cataract Surgery of Intraocular Lens Constant Personalization using the Haigis Formula." *Journal of Cataract and Refractive Surgery* 36, no. 7 (2010): 1081-1089. doi:10.1016/j.jcrs.2009.12.050.

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
**Description**

Model	ANKORIS						
Material	26% Hydrophilic Acrylic						
Overall diameter	11.40mm						
Optic diameter	6.00mm						
Optic	Biconvex Aspheric Monofocal						
Haptic design	Double C-loop & Posterior Angulated Haptic						
Filtration	UV & Blue Light						
Refractive index	1.46						
Abbe number	58						
Injection system	Medical Accuject 2.0 up to 24.5D & Medical Accuject 2.1/2.2 up to 30D						
Spherical power <sup>1</sup>	+6D to +30D (0.5D steps)						
Cylinder power (IOL plane) <sup>2</sup>	1.50 - 2.25 - 3.00 - 3.75 - 4.50 - 5.25 - 6.00						
Suggested A constant <sup>3</sup>	Interferometry						
	Hoffer Q: pACD		5.59				
	Holladay 1: Sf		1.83				
	Barrett: LF		1.86				
	SRK/T: A		118.95				
	Haigis <sup>4</sup> : a0; a1; a2		1.36; 0.4; 0.1				
Cylinder power at IOL plane	ANKORIS 1.5	ANKORIS 2.25	ANKORIS 3.0	ANKORIS 3.75	ANKORIS 4.5	ANKORIS 5.25	ANKORIS 6.0
	1.50D	2.25D	3.00D	3.75D	4.50D	5.25D	6.00D
Cylinder power at corneal plane <sup>4</sup>	1.03D	1.55D	2.06D	2.57D	3.08D	3.60D	4.11D

<sup>1</sup> Values estimated only; surgeons are recommended to personalize their A-constant based on their surgical techniques and equipment, experience with the lens model and postoperative results.

<sup>2</sup> Not optimized.

<sup>3</sup> Savini G., J Cataract Refract Surg 2013; 39:1900-1903.

<sup>4</sup> Please check the availability of spherical and cylinder powers with your sales representative.

Note: The ANKORIS intraocular lens is not FDA approved.

Contact information:  
[www.bvimedical.com/customer-support/](http://www.bvimedical.com/customer-support/)

[bvimedical.com](http://bvimedical.com)

BVI and all other trademarks (unless noted otherwise) are property of BVI. BVI ©2022

590626-03

## Product Information

Manufacturer	PhysIOL s.a. - Liège Science Park Allée des Noisetiers 4 B-4031 Belgium +32 4 361 05 49 <a href="mailto:physiol@bvimedical.com">physiol@bvimedical.com</a>
Certificate information	CE (EU) 2017/745, Annex IX Chapter II : MDR 735733 R000 QMS (EU) 2017/745, Annex IX Chapter I and III : MDR 735719 R000 ISO 13485:2016 & EN ISO 13485:2016 : MD 658518 ISO 13485:2016 : MDSAP 691544
Shelf life	Five (5) years from manufacturing date
Intended Purpose	The posterior chamber intraocular lens which is intended to be placed into the capsular bag for the replacement of the human lens to achieve the visual correction of aphakia in adult patients in whom the cataractous lens has been removed.
Indication for use	The lens should be used as intended in adult patients, with pre-existing astigmatism, surgically treated for cataract, who desire improved uncorrected far vision, with reduced spectacle dependence.
Product Composition	No products of animal or human origin are present in the implant. The intraocular lens is 100% composed of the covalently crosslinked medical quality material HELIOFLEX, which is a (2-hydroxyethylmethacrylate; methylmethacrylate) copolymer including a UV and a blue light-filtering chromophores covalently bound to the material.
For sterile product	All IOLs from PhysIOL are steam sterilized
Packaging Material	Holder (Polypropylene) Container (Polypropylene) Storage liquid (0.9% NaCl solution) Aluminium lid (Aluminium Gold) Container label (paper) Blister PP (Polypropylene) Tyvek lid
Product Class	Classified as Class IIb implantable long-term surgically invasive medical devices under Rule 8 of Annex VIII of MDR 2017/745. Not available in the United States

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# B! | FLEX T

677TA | 677TAY

Posterior-chamber toric intraocular lens, pseudophakic, hydrophilic, with optional blue light filter

 CLEAR IOL

 YELLOW IOL

POWER RANGE				
Standard Powers	0.0 D + +30.0 D (0.5 D steps)			
Extreme Powers	-10.0 D + -1.0 D (1.0 D steps) +31.0 D + +35.0 D (1.0 D steps)			
Cylinders*	1.0 D	1.5 D	2.25 D	3.0 D
	3.75 D	4.5 D	5.25 D	6.0 D
MATERIAL				
Copolymer of hydrophobic and hydrophilic monomers with 25% water content, UV blocker and optional blue light filter				
Refractive Index	1.46			
Abbe Number	58			
DESIGN				
Optic Diameter Overall Diameter	6 mm		13 mm	
Optic Design	Aspheric - Aberration Neutral Convex-concave (-10.0 D + +5.0 D) Biconvex (+5.5 D + +35.0 D)			
Nominal A-constant	118.9			
Haptic Design	Posterior vaulting fenestrated C-loops with 0° angulation			
PCO Prevention	360° sharp edge			
Shelf Life	5 years after sterilization			
INJECTION				
Estimated Incision Size	2.2 - 2.4 mm			
Recommended injector	Compatibility chart available: <a href="https://www.medicontur.com/professionals/compatibility">medicontur.com/professionals/compatibility</a>			
CONSTANTS <i>for optical and immersion ultrasound measurements</i>				
SRK/T A-Constant	Haigis		Hoffer Q	
A = 118.828** (optimized)	a0 = 0.190** a1 = 0.192** a2 = 0.173**		pACD = 5.431**	
Holladay I	Barrett Universal II		Holladay II	
SF = 1.682**	LF = 1.79		ACD = 5.45 SF = 1.69	

\*Extreme cylinders available: +6.75 D, +7.5 D, +8.25 D, +9.0 D, +10.0 D  
\*\* Optimized constants based on more than 400 patients

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## Apendice 6. LIO Medicontur ELON

877PEY

Single piece, sterile, foldable aspheric hydrophobic acrylic intraocular lens (IOL) with UV blocker and blue light filter. The optic is EDOF. The IOL is available as part of a preloaded injector system.

YELLOW IOL


The safety and performance of Medicontur IOLs have been tested in clinical investigations and the clinical data demonstrate that these IOLs meet the requirements of the state of the art.

**MEDICONTUR**

POWER RANGE		
Standard Powers	+8.0 D + +30.0 D (0.5 D steps)	
Extreme Powers	+31.0 D + +35.0 D (1.0 D steps)	
MATERIAL		
Hydrophobic acrylic with UV blocker and blue light filter		
Refractive Index	1.47	
Abbe Number	58	
DESIGN		
Optic Diameter Overall Diameter	6 mm	13 mm
Optic Design	Aspheric - Aberration Neutral Biconvex (+8.0 D + +35.0 D)	
Nominal A-constant	118.9	
EDOF Zone	Refraction-based central Wavefront Linking zones	
Haptic Design	Posterior vaulting fenestrated C-loops with 0° angulation	
PCO Prevention	360° sharp edge	
Shelf Life	2.5 years after sterilization	
INJECTION		
Estimated Incision Size	2.2 mm	
Compatible injector	POB-MA SINGLE-USE PRELOADED INJECTOR	
CONSTANTS		
for optical and immersion ultrasound measurements		
SRK/T A-Constant	Haigis	Hoffer Q
A = 118.9	a0 = 1.32 a1 = 0.4 a2 = 0.1	pACD = 5.46
Holladay I	Barrett Universal II	Holladay II
SF = 1.7	LF = 1.83	ACD = 5.49 SF = 1.73

IMA-IOL-ELONG-IMD02S-DTS-EN 202008

Apéndice 7: Documentación LIO Liberty (IOL Liberty 677MY. Medicontur Medical Engineering Ltd., Zsámbék, Hungary)

## LIBERTY 677MY

Single piece, sterile, foldable acrylic aspheric hydrophilic posterior chamber intraocular lens (IOL) with UV blocker and blue light filter. The IOL is non-preloaded, which has to be manually loaded into a compatible injector.

<b>Powers:</b>	+8.0 D → +35.0 D (0.5 D steps)
<b>Addition (at the IOL plane):</b>	+3.5 D near +1.75 D intermediate
<b>Material:</b>	Copolymer of hydrophobic and hydrophilic monomers with 25% water content, UV blocker and blue light filter
<b>RI:</b>	1.46
<b>Abbe No:</b>	58
<b>Optic Diameter:</b>	6 mm
<b>Overall Diameter:</b>	13 mm
<b>Optic Design:</b>	Aspheric - Aberration Neutral Biconvex (+8.0 D → +35.0 D)
<b>Nominal A-constant:</b>	118.9
<b>Diffractive Zone:</b>	3 mm diameter diffractive array on the anterior surface utilizing EPS technology
<b>Haptic Design:</b>	Posterior vaulting fenestrated C-loops with 0° angulation
<b>PCO Prevention:</b>	360° sharp edge
<b>Estimated Incision Size:</b>	1.8 - 2.2 mm
<b>Shelf Life:</b>	5 years after sterilization
<b>Recommended Injector:</b>	Compatibility chart available: <a href="https://medicontur.com/professionals/compatibility">https://medicontur.com/professionals/compatibility</a>

Global site

EC Certificate Full Quality Assurance System:  
Certificate HU19/8471

The management system of

## **MEDICONTUR Medical Engineering Ltd.**

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has been assessed and certified as meeting the requirements of

**Directive 93/42/EEC**  
on medical devices, Annex II (excluding Section 4)

For the following products

**The scope of registration appears on page 2 of this certificate**

This certificate is valid from 27 April 2021 until 02 May 2023  
and remains valid subject to satisfactory surveillance audits.  
Issue 3. Certified since 20 May 1997.

Certification is based on reports numbered HU/BUD HU0164QYH

This is a multi-site certification.  
Additional site details are listed on the subsequent page.

Authorised by

Global Medical Devices Head of Notified Body

**SGS Belgium NV, Notified Body 1639**

SGS House Noorderlaan 87 2030 Antwerp Belgium  
t +32 (0)3 545-48-48 f +32 (0)3 545-48-49 [www.sgs.com](http://www.sgs.com)

LPM05007 - Certificate CE1639 Annex II-4\_EN rev. 02

Page 1 of 2

Certificate HU19/8471 continued

# **MEDICONTUR Medical Engineering Ltd.**

## **Directive 93/42/EEC**

on medical devices, Annex II (excluding Section 4)  
Issue 3

Detailed scope

**Sterile single use foldable lens injectors**  
**Disposable Injection Kits**  
**Sterile Viscoelastic solutions for ophthalmic use**  
**Sterile Capsular Tension Ring Preloaded In A Single Use Injection Kit**  
**Aspheric hydrophilic acrylic IOL for implantation to the capsular bag**  
**Aspheric hydrophilic acrylic IOL with blue light filter**  
**for implantation to the capsular bag**  
**Hydrophilic acrylic IOL for implantation into the ciliary sulcus**  
**Aspheric hydrophobic acrylic IOL for implantation into the capsular bag**  
**Aspheric hydrophobic acrylic IOL with blue light filter**  
**for implantation into the capsular bag**

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Where the above scope includes class III medical device(s), a valid EC Design Examination Certificate according to Annex II (Section 4) is a mandatory requirement for each device in addition to this certificate to place that device on the market.

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Trifocal  
Hydrophilic

### Technical Specifications

Commercial name	MICRO F		POD F	
Material	25% hydrophilic acrylic		26% hydrophilic acrylic	
Overall diameter	10.75 mm		11.40 mm	
Optic diameter	6.15 mm		6.00 mm	
Optic	Biconvex aspheric (-0.11μ SA) trifocal diffractive FINEVISION			
Filtration	UV & blue light			
Refractive index	1.46			
Abbe number	58			
Angulation	5°			
Additional power	+ 1.75D for intermediate vision and + 3.50D for near vision			
Injection system	Medicel Viscoject Bio 1.8/2.2 up to 24.5D Medicel Accuject 1.8/2.0 up to 24.5D Medicel Accuject 2.1/2.2 up to 35D		Medicel Accuject 2.0 up to 24.5D Medicel Accuject 2.1/2.2 up to 35D	
Incision size	≥ 1.8 mm		≥ 2.0 mm	
Spherical power	10D to 35D (0.5D steps)		6D to 35D (0.5D steps)	
Square edge	360°			
Nominal manufacturer A constant	118.80		118.95	
Suggested A constant <sup>1</sup>	Interferometry	Ultrasound	Interferometry	Ultrasound
Hoffer Q: pACD	5.35	5.26	5.59	5.35
Holladay 1: Sf	1.60	1.48	1.83	1.57
Barrett: LF	1.78	-	1.86	-
SRK/T: A	118.80	118.59	118.95	118.73
Haigis <sup>2</sup> : a0; a1; a2	1.36; 0.4; 0.1	1.13; 0.4; 0.1	1.36; 0.4; 0.1	1.13; 0.4; 0.1

<sup>1</sup> Estimates only; surgeons are recommended to use their own values based upon their personal experience. Refer to our website for updates.

<sup>2</sup> Not optimized.

**Contact Information:**  
[www.bvimedical.com/customer-support/](http://www.bvimedical.com/customer-support/)

[bvimedical.com](http://bvimedical.com)

BVI and all other trademarks (unless noted otherwise) are property of BVI. BVI ©2020

590627-01

## Product Information

	MICRO F	POD F
Manufacturer	PhysIOL s.a. - Liège Science Park Allée des Noisetiers 4 B-4031 Belgium +32 4 361 05 49 physiol@bvmedical.com	PhysIOL s.a. - Liège Science Park Allée des Noisetiers 4 B-4031 Belgium +32 4 361 05 49 physiol@bvmedical.com
Certificate information	CE: Certificate N° CE658516 ISO 13485:2016: Certificate n° MD658518 MDSAP: Certificate N° MDSAP 691544 ISO 9001:2015: Certificate N° FM 658519	CE: Certificate N° CE658516 ISO 13485:2016: Certificate n° MD658518 MDSAP: Certificate N° MDSAP 691544 ISO 9001:2015: Certificate N° FM 658519
Shelf life	Five (5) years from manufacturing date	Five (5) years from manufacturing date
Intended Use	Intended use (for all IOLs): The posterior chamber intraocular lens which is intended to be placed into the capsular bag for the replacement of the human lens to achieve the visual correction of aphakia in adult patients in whom the cataractous lens has been removed by extracapsular cataract extraction.	Intended use (for all IOLs): The posterior chamber intraocular lens which is intended to be placed into the capsular bag for the replacement of the human lens to achieve the visual correction of aphakia in adult patients in whom the cataractous lens has been removed by extracapsular cataract extraction.
Indication for use	The lens should be used as intended in patients surgically treated for cataract, with possibly associated presbyopia, who desire improved uncorrected far vision, useful near and intermediate visual functions and reduced spectacle dependence.	The lens should be used as intended in patients surgically treated for cataract, with possibly associated presbyopia, who desire improved uncorrected far vision, useful near and intermediate visual functions and reduced spectacle dependence.
Product Composition	No products of animal or human origin are present in the implant. The implant is made of the HELIO25 material, composed of an acrylate copolymer Hydroxyethyl methacrylate (HEMA) and Ethoxyethyl methacrylate (EOEMA), including a UV and blue light filter	No products of animal or human origin are present in the implant. The implant is made of the HELIOFLEX material, composed of an acrylate copolymer Hydroxyethyl methacrylate (HEMA) and Methyl methacrylate (MMA), including a UV and blue light filter
For sterile product	All IOLs from PhysIOL are steam sterilized	All IOLs from PhysIOL are steam sterilized
Packaging Material	Holder (Polypropylene) Container (Polypropylene) Storage liquid (0.9% NaCl solution) Aluminium lid (Aluminium Gold) Container label (paper) Blister PP (Polypropylene) Tyvek lid	Holder (Polypropylene) Container (Polypropylene) Storage liquid (0.9% NaCl solution) Aluminium lid (Aluminium Gold) Container label (paper) Blister PP (Polypropylene) Tyvek lid
Product Class	MDD Class IIb Sterile, According to European Medical Device Directive 93/42/EEC Not available in the United States	MDD Class IIb Sterile, According to European Medical Device Directive 93/42/EEC Not available in the United States

Contact Information:  
[www.bvmedical.com/customer-support/](http://www.bvmedical.com/customer-support/)

BVI and all other trademarks (unless noted otherwise) are property of BVI. BVI ©2020

590627-01

## Injection Guidelines

The Medical Viscoject Bio 1.8/2.2 and Accuject 1.8/2.0/2.1/2.2 injection system are recommended for implanting the FINEVISION lenses.

This fully single-use system represents total reliability for safe and effective lens injections.


Its compact design with integrated cartridge enables a simple, predictable loading and positioning of the lens.

### FINEVISION MICRO F

Medical Viscoject Bio 1.8/2.2 up to 24.5D

Medical Accuject 1.8/2.0 up to 24.5D and 2.1/2.2 up to 35D

Guidelines with Viscoject:


1. Apply ophthalmic viscoelastic device (OVD) into the tip and the loading chamber of the injector cartridge.
2. Remove the lens from the lens holder. Position the lens into the cartridge in such way that the two haptics with the notches are pointing at 1 and 7 o'clock.
3. Exert slight pressure onto the lens optic and make sure that all haptics are inside before further closing the cartridge.
4. Close the cartridge and check the position of the lens. Once the "click-lock" mechanism engages, the lens is securely loaded.
5. Engage the cartridge in the injector.
6. Press the injector plunger forward and push the lens into the conical tip of the cartridge. Pull the plunger back a few millimeters and then inject the lens in one continuous motion. For gentle implantation, it is not necessary to fully push the plunger to the bottom of the cartridge.

### FINEVISION POD F

Accuject 2.0 for lens diopters up to 24.5D

Accuject 2.1/2.2 for lens diopters up to 35D

Guidelines with Accuject:

1. Apply ophthalmic viscoelastic device (OVD) into the tip and the loading chamber of the injector cartridge.
2. Remove the lens from the lens holder. Position the lens into the cartridge in such way that the two haptics with the notches are pointing at 1 and 7 o'clock.
3. Exert slight pressure onto the lens optic and make sure that all haptics are inside before further closing the cartridge. Close the cartridge and check the position of the lens.
4. Once the "click-lock" mechanism engages, the lens is securely loaded and ready for injection.
5. Press the injector plunger forward and push the lens into the conical tip of the cartridge.
6. Pull the plunger back a few millimeters and then inject the lens in one continuous motion. For gentle implantation, it is not necessary to fully push the plunger to the bottom of the cartridge.

Contact Information:  
[www.bvimedical.com/customer-support/](http://www.bvimedical.com/customer-support/)

BVI and all other trademarks (unless noted otherwise) are property of BVI. BVI ©2020

590627-01